



SOFTWARE RELEASE NOTES

MCS-CONNECT Version 18.00B 02/16/2016

Updates from previous versions:

New Features:

1. Add the new folders for custom_algo_chiller.c, custom_algo_chiller.h, custom_algo_rtu.c and custom_algo_rtu.h.
2. Add protocol option for the rs485 #2 for future use.
3. Add two new Setpoint STPT_COMPRESSOR_5_FLA (# 86) and STPT_COMPRESSOR_6_FLA (#87) for the added compressors.
4. Add new element to MACHINE_EEPROM_TYEPDEF structure Bytes # 7-8 for Bacnet Device ID.

Bug Fixes:

1. When checking for MOP the logic was entering the EXV operating state when in Pre Pump Down, setting the EXV to the 10%. By checking for the state to not be in pre pump down before enter the MOP logic solve this bug.
2. The evap. fan control reverse target is now controled from config. using Control Method bit.
3. Logic for reverse valve was resetting the Compressor RO. The reverse valve logic is only use with SKM SW, by adding statement to use only with SKM SW solve the bug.

Any questions regarding this release, contact:
support@mcscontrols.com

Micro Control Systems, Inc. 5580 Enterprise Parkway Fort Myers, Florida 33905
(239)694-0089 FAX: (239)694-0031 www.mcscontrols.com

Information contained in this document has been prepared by Micro Control Systems, Inc. and is copyright © protected 2018. Copying or distributing this document is prohibited unless expressly approved by MCS.